

Product datasheet for **TP720336L**

SLITL2 (VASN) (NM_138440) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human vasorin (VASN) |
| Species: | Human |
| Expression Host: | HEK293 |
| Expression cDNA Clone or AA Sequence: | Cys24-Pro575 |
| Tag: | C-His |
| Predicted MW: | 60.7 kDa |
| Concentration: | lot specific |
| Purity: | >95% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl |
| Endotoxin: | < 0.1 EU per µg protein as determined by LAL test |
| Reconstitution Method: | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |
| Storage: | Store at -80°C. |
| Stability: | Stable for at least 6 months from date of receipt under proper storage and handling conditions. |
| RefSeq: | NP_612449 |
| Locus ID: | 114990 |
| UniProt ID: | Q6EMK4 |
| Cytogenetics: | 16p13.3 |
| Synonyms: | SLITL2 |
| Summary: | May act as an inhibitor of TGF-beta signaling.[UniProtKB/Swiss-Prot Function] |
| Protein Families: | Secreted Protein, Transmembrane |



[View online »](#)

Product images:

