

Product datasheet for **TP720347L**

CD93 (NM_012072) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CD93 molecule (CD93)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Ala24-Lys580
Tag:	C-His
Predicted MW:	59.2 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_036204
Locus ID:	22918
UniProt ID:	Q9NPY3
Cytogenetics:	20p11.21
Synonyms:	C1qR(P); C1QR1; C1qRP; CDw93; dj737E23.1; ECSM3; MXRA4



[View online »](#)

Summary:

The protein encoded by this gene is a cell-surface glycoprotein and type I membrane protein that was originally identified as a myeloid cell-specific marker. The encoded protein was once thought to be a receptor for C1q, but now is thought to instead be involved in intercellular adhesion and in the clearance of apoptotic cells. The intracellular cytoplasmic tail of this protein has been found to interact with moesin, a protein known to play a role in linking transmembrane proteins to the cytoskeleton and in the remodelling of the cytoskeleton. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome, ES Cell Differentiation/IPS, Transmembrane

Product images: