

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TP720357

### RAGE (AGER) (NM\_001136) Human Recombinant Protein

# **Product data:**

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human advanced glycosylation end product-specific receptor (AGER), transcript variant 1
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Ala23-Ala344
Tag:	C-His
Predicted MW:	35.2 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 $\mu m$ filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Bioactivity:	Binding assay (PMID: <u>29420566</u> )
Endotoxin:	< 0.1 EU per $\mu$ g protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<u>NP 001127</u>
Locus ID:	177
UniProt ID:	<u>Q15109, A0A1U9X785, B4DNX3</u>
Cytogenetics:	6p21.32
Synonyms:	RAGE; SCARJ1; sRAGE



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

# SAGE (AGER) (NM\_001136) Human Recombinant Protein – TP720357

Summary: The advanced glycosylation end product (AGE) receptor encoded by this gene is a member of the immunoglobulin superfamily of cell surface receptors. It is a multiligand receptor, and besides AGE, interacts with other molecules implicated in homeostasis, development, and inflammation, and certain diseases, such as diabetes and Alzheimer's disease. Many alternatively spliced transcript variants encoding different isoforms, as well as non-proteincoding variants, have been described for this gene (PMID:18089847). [provided by RefSeq, May 2011]

#### Protein Families: Druggable Genome, Secreted Protein, Transmembrane

#### **Product images:**



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US