

## Product datasheet for **TP720370XL**

### Cathepsin Z (CTSZ) (NM\_001336) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cathepsin Z (CTSZ)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Gly24-Val303
Tag:	C-His
Predicted MW:	32.5 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001327</a>
Locus ID:	1522
UniProt ID:	<a href="#">Q9UBR2</a>
Cytogenetics:	20q13.32
Synonyms:	CTSX
Summary:	The protein encoded by this gene is a lysosomal cysteine proteinase and member of the peptidase C1 family. It exhibits both carboxy-monopeptidase and carboxy-dipeptidase activities. The encoded protein has also been known as cathepsin X and cathepsin P. This gene is expressed ubiquitously in cancer cell lines and primary tumors and, like other members of this family, may be involved in tumorigenesis. [provided by RefSeq, Oct 2008]
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Lysosome



[View online »](#)

**Product images:**

