

Product datasheet for **TP720371M**

Nectin 2 (NECTIN2) (NM_001042724) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human poliovirus receptor-related 2 (herpesvirus entry mediator B) (PVRL2), transcript variant delta
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Gln32-Leu360
Tag:	C-His
Predicted MW:	36.6 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_001036189
Locus ID:	5819
UniProt ID:	Q92692
Cytogenetics:	19q13.32
Synonyms:	CD112; HVEB; PRR2; PVRL2; PVRR2



[View online »](#)

Summary:

This gene encodes a single-pass type I membrane glycoprotein with two Ig-like C2-type domains and an Ig-like V-type domain. This protein is one of the plasma membrane components of adherens junctions. It also serves as an entry for certain mutant strains of herpes simplex virus and pseudorabies virus, and it is involved in cell to cell spreading of these viruses. Variations in this gene have been associated with differences in the severity of multiple sclerosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome, Transmembrane

Protein Pathways:

Adherens junction, Cell adhesion molecules (CAMs)

Product images: