

## Product datasheet for **TP720391L**

### Ephrin B2 (EFNB2) (NM\_004093) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ephrin-B2 (EFNB2)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Ile28-Ala229
Tag:	C-His
Predicted MW:	23.2 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_004084</a>
Locus ID:	1948
UniProt ID:	<a href="#">P52799</a>
Cytogenetics:	13q33.3
Synonyms:	EPLG5; Htk-L; HTKL; LERK5



[View online »](#)

**Summary:**

This gene encodes a member of the ephrin (EPH) family. The ephrins and EPH-related receptors comprise the largest subfamily of receptor protein-tyrosine kinases and have been implicated in mediating developmental events, especially in the nervous system and in erythropoiesis. Based on their structures and sequence relationships, ephrins are divided into the ephrin-A (EFNA) class, which are anchored to the membrane by a glycosylphosphatidylinositol linkage, and the ephrin-B (EFNB) class, which are transmembrane proteins. This gene encodes an EFNB class ephrin which binds to the EPHB4 and EPHA3 receptors. [provided by RefSeq, Jul 2008]

**Protein Families:**

Druggable Genome, Transmembrane

**Protein Pathways:**

Axon guidance

**Product images:**