

## Product datasheet for **TP720416XL**

### CDw75 (ST6GAL1) (NM\_003032) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ST6 beta-galactosamide alpha-2,6-sialyltransferase 1 (ST6GAL1), transcript variant 2
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Lys27-Cys406
Tag:	C-His
Predicted MW:	44.6 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_003023</a>
Locus ID:	6480
UniProt ID:	<a href="#">P15907</a> , <a href="#">B2R513</a>
Cytogenetics:	3q27.3
Synonyms:	SIAT1; ST6Gall; ST6N
Summary:	This gene encodes a member of glycosyltransferase family 29. The encoded protein is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The protein, which is normally found in the Golgi but can be proteolytically processed to a soluble form, is involved in the generation of the cell-surface carbohydrate determinants and differentiation antigens HB-6, CD75, and CD76. This gene has been incorrectly referred to as CD75. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2017]



[View online »](#)

Protein Families: Secreted Protein

Protein Pathways: Metabolic pathways, N-Glycan biosynthesis

**Product images:**

