

Product datasheet for **TP720435L**

PDILT (NM_174924) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein disulfide isomerase-like, testis expressed (PDILT)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Ser21-Leu584
Tag:	C-His
Predicted MW:	65.5 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_777584
Locus ID:	204474
UniProt ID:	Q8N807
Cytogenetics:	16p12.3
Synonyms:	PDIA7
Summary:	This gene encodes a member of the disulfide isomerase (PDI) family of endoplasmic reticulum (ER) proteins that catalyze protein folding and thiol-disulfide interchange reactions. The encoded protein has an N-terminal ER-signal sequence, two thioredoxin (TRX) domains with non-classical Ser-Lys-Gln-Ser and Ser-Lys-Lys-Cys motifs, respectively, two TRX-like domains, and a C-terminal ER-retention sequence. The protein lacks oxidoreductase activity in vitro and probably functions as a chaperone. This gene's expression appears to be limited to the testis. [provided by RefSeq, Dec 2016]



[View online »](#)

Protein Families: Druggable Genome

Product images:

