

Product datasheet for **TP720491**

Nogo Receptor (RTN4R) (NM_023004) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human reticulon 4 receptor (RTN4R)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Cys27-Ser447
Tag:	C-His
Predicted MW:	46.3 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_075380
Locus ID:	65078
UniProt ID:	Q9BZR6
Cytogenetics:	22q11.21
Synonyms:	NGR; NOGOR
Summary:	This gene encodes the receptor for reticulon 4, oligodendrocyte myelin glycoprotein and myelin-associated glycoprotein. This receptor mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system. [provided by RefSeq, Jul 2008]



[View online »](#)

Protein Families: Druggable Genome

Product images:

