

Product datasheet for TP720492XL

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Semaphorin 3c (SEMA3C) (NM_006379) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human sema domain, immunoglobulin domain (Ig), short basic

domain, secreted, (semaphorin) 3C (SEMA3C)

Species: Human Expression Host: HEK293

Expression cDNA Clone

or AA Sequence:

Gly21-Ser738

Tag: C-His

Predicted MW: 82.4 kDa

Concentration: lot specific

Purity: >95% as determined by reducing SDS-PAGE

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: < 0.1 EU per µg protein as determined by LAL test

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 006370

 Locus ID:
 10512

 UniProt ID:
 Q99985

 Cytogenetics:
 7q21.11

Synonyms: SEMAE; SemE





Summary:

This gene encodes a secreted glycoprotein that belongs to the semaphorin class 3 family of neuronal guidance cues. The encoded protein contains an N-terminal sema domain, integrin and immunoglobulin-like domains, and a C-terminal basic domain. Homodimerization and proteolytic cleavage of the C-terminal propeptide are necessary for the function of the encoded protein. It binds a neuropilin co-receptor before forming a heterotrimeric complex with an associated plexin. An increase in the expression of this gene correlates with an increase in cancer cell invasion and adhesion. Naturally occurring mutations in this gene are associated with Hirschsprung disease. [provided by RefSeq, May 2017]

Protein Families: Secreted Protein
Protein Pathways: Axon guidance

Product images:

