

## Product datasheet for **TP720495XL**

### VMO1 (NM\_001144940) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human vitelline membrane outer layer 1 homolog (chicken) (VMO1), transcript variant 3
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Gln25-Ser202
Tag:	C-His
Predicted MW:	20.1 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001138412</a>
Locus ID:	284013
UniProt ID:	<a href="#">Q7Z5L0</a>
Cytogenetics:	17p13.2
Synonyms:	ERGA6350; PRO21055
Protein Families:	Secreted Protein



[View online »](#)

**Product images:**

