

Product datasheet for **TP720507M**

PBLD (NM_022129) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phenazine biosynthesis-like protein domain containing (PBLD), transcript variant 1
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Ala288
Tag:	N-His
Predicted MW:	33.9 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 30% Glycerol, pH 8.0.
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_071412
Locus ID:	64081
UniProt ID:	P30039 , A0A024QZK5
Cytogenetics:	10q21.3
Synonyms:	HEL-S-306; MAWBP; MAWDBP



[View online »](#)

Product images:

