

## Product datasheet for **TP720507XL**

### **PBLD (NM\_022129) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant protein of human phenazine biosynthesis-like protein domain containing (PBLD), transcript variant 1
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	Met1-Ala288
<b>Tag:</b>	N-His
<b>Predicted MW:</b>	33.9 kDa
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>95% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 30% Glycerol, pH 8.0.
<b>Endotoxin:</b>	< 0.1 EU per µg protein as determined by LAL test
<b>Storage:</b>	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Stability:</b>	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
<b>RefSeq:</b>	<a href="#">NP_071412</a>
<b>Locus ID:</b>	64081
<b>UniProt ID:</b>	<a href="#">P30039</a> , <a href="#">A0A024QZK5</a>
<b>Cytogenetics:</b>	10q21.3
<b>Synonyms:</b>	HEL-S-306; MAWBP; MAWDBP



[View online »](#)

**Product images:**

