

Product datasheet for **TP720509XL**

BDH2 (NM_020139) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human 3-hydroxybutyrate dehydrogenase, type 2 (BDH2)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Leu245
Tag:	N-His
Predicted MW:	28.9 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_064524
Locus ID:	56898
UniProt ID:	Q9BUT1
Cytogenetics:	4q24
Synonyms:	DHRS6; EFA6R; PRO20933; SDR15C1; UCPA-OR; UNQ6308
Summary:	Dehydrogenase that mediates the formation of 2,5-dihydroxybenzoic acid (2,5-DHBA), a siderophore that shares structural similarities with bacterial enterobactin and associates with LCN2, thereby playing a key role in iron assimilation and homeostasis. Plays a role in susceptibility to bacterial infection by providing an assimilable source of iron that is exploited by pathogenic bacteria (By similarity). Also acts as a 3-hydroxybutyrate dehydrogenase (PubMed:16380372).[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Butanoate metabolism, Metabolic pathways, Synthesis and degradation of ketone bodies

Product images: