

Product datasheet for **TP720513L**

PITPN (PITPNA) (NM_006224) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphatidylinositol transfer protein, alpha (PITPNA)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Asp270
Tag:	N-His
Predicted MW:	34 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 1mM EDTA, 1mM DTT, pH 8.0.
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_006215
Locus ID:	5306
UniProt ID:	Q00169 , V9HWC5
Cytogenetics:	17p13.3
Synonyms:	HEL-S-36; PI-TPalpha; PITPN; VIB1A
Summary:	This gene encodes a member of a family of lipid-binding proteins that transfer molecules of phosphatidylinositol or phosphatidylcholine between membrane surfaces. The protein is implicated in phospholipase C signaling and in the production of phosphatidylinositol 3,4,5-trisphosphate (PIP3) by phosphoinositide-3-kinase.[provided by RefSeq, Sep 2009]



[View online »](#)

Product images:

