

Product datasheet for **TP720517L**

CLIC4 (NM_013943) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chloride intracellular channel 4 (CLIC4), nuclear gene encoding mitochondrial protein
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Lys253
Tag:	N-His
Predicted MW:	30.9 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_039234
Locus ID:	25932
UniProt ID:	Q9Y696 , Q6FIC5
Cytogenetics:	1p36.11
Synonyms:	CLIC4L; H1; huH1; MTCLIC; p64H1
Summary:	Chloride channels are a diverse group of proteins that regulate fundamental cellular processes including stabilization of cell membrane potential, transepithelial transport, maintenance of intracellular pH, and regulation of cell volume. Chloride intracellular channel 4 (CLIC4) protein, encoded by the CLIC4 gene, is a member of the p64 family; the gene is expressed in many tissues and exhibits a intracellular vesicular pattern in Panc-1 cells (pancreatic cancer cells). [provided by RefSeq, Jul 2008]



[View online »](#)

Protein Families: Druggable Genome, Ion Channels: Other

Product images:

