

## **Product datasheet for TP720517XL**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CLIC4 (NM\_013943) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human chloride intracellular channel 4 (CLIC4), nuclear gene

encoding mitochondrial protein

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

Met1-Lys253

Tag: N-His

**Predicted MW:** 30.9 kDa **Concentration:** lot specific

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

**Endotoxin:** < 0.1 EU per µg protein as determined by LAL test

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

**RefSeq:** NP 039234

**Locus ID:** 25932

 UniProt ID:
 Q9Y696, Q6FIC5

Cytogenetics: 1p36.11

Synonyms: CLIC4L; H1; huH1; MTCLIC; p64H1

**Summary:** Chloride channels are a diverse group of proteins that regulate fundamental cellular

processes including stabilization of cell membrane potential, transepithelial transport, maintenance of intracellular pH, and regulation of cell volume. Chloride intracellular channel 4 (CLIC4) protein, encoded by the CLIC4 gene, is a member of the p64 family; the gene is

expressed in many tissues and exhibits a intracellular vesicular pattern in Panc-1 cells

(pancreatic cancer cells). [provided by RefSeg, Jul 2008]





**Protein Families:** Druggable Genome, Ion Channels: Other

## **Product images:**

