

## Product datasheet for **TP720525XL**

### Inosine triphosphate pyrophosphatase (ITPA) (NM\_033453) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human inosine triphosphatase (nucleoside triphosphate pyrophosphatase) (ITPA), transcript variant 1
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ala2-Ala194
Tag:	C-His
Predicted MW:	22.5 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_258412</a>
Locus ID:	3704
UniProt ID:	<a href="#">Q9BY32</a> , <a href="#">A0A0S2Z3W7</a>
Cytogenetics:	20p13
Synonyms:	C20orf37; DEE35; dj79416.3; HLC14-06-P; ITPase; My049; NTPase
Summary:	This gene encodes an inosine triphosphate pyrophosphohydrolase. The encoded protein hydrolyzes inosine triphosphate and deoxyinosine triphosphate to the monophosphate nucleotide and diphosphate. This protein, which is a member of the HAM1 NTPase protein family, is found in the cytoplasm and acts as a homodimer. Defects in the encoded protein can result in inosine triphosphate pyrophosphorylase deficiency which causes an accumulation of ITP in red blood cells. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jun 2012]



[View online »](#)

**Protein Families:** Druggable Genome

**Protein Pathways:** Drug metabolism - other enzymes, Metabolic pathways, Purine metabolism, Pyrimidine metabolism

**Product images:**

