

Product datasheet for TP720527XL

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GADD45G (NM_006705) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human growth arrest and DNA-damage-inducible, gamma

(GADD45G)

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

Met1-Glu159

Tag: N-His

Predicted MW: 19.3 kDa

Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: < 0.1 EU per µg protein as determined by LAL test

Reconstitution Method: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

Stability: Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 006696

 Locus ID:
 10912

 UniProt ID:
 095257

 Cytogenetics:
 9q22.2

Synonyms: CR6; DDIT2; GADD45gamma; GRP17





Summary:

This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. The protein encoded by this gene responds to environmental stresses by mediating activation of the p38/JNK pathway via MTK1/MEKK4 kinase. The GADD45G is highly expressed in placenta. [provided by RefSeq, Jul 2008]

Protein Pathways:

Cell cycle, MAPK signaling pathway, p53 signaling pathway

Product images:

