

## Product datasheet for **TP720551L**

### ZNF75A (NM\_153028) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger protein 75a (ZNF75A)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	S58-K162
Tag:	N-6His
Predicted MW:	14.3 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_694573</a>
Locus ID:	7627
UniProt ID:	<a href="#">Q96N20</a> , <a href="#">Q68CU0</a>
Cytogenetics:	16p13.3
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors



[View online »](#)

**Product images:**

