

Product datasheet for **TP720556XL**

Glutathione Synthetase (GSS) (NM_000178) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glutathione synthetase (GSS)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ala2-Val474
Tag:	C-His
Predicted MW:	53.5 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_000169
Locus ID:	2937
UniProt ID:	P48637 , V9HWJ1
Cytogenetics:	20q11.22
Synonyms:	GSHS; HEL-S-64p; HEL-S-88n
Summary:	Glutathione is important for a variety of biological functions, including protection of cells from oxidative damage by free radicals, detoxification of xenobiotics, and membrane transport. The protein encoded by this gene functions as a homodimer to catalyze the second step of glutathione biosynthesis, which is the ATP-dependent conversion of gamma-L-glutamyl-L-cysteine to glutathione. Defects in this gene are a cause of glutathione synthetase deficiency. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Glutathione metabolism, Metabolic pathways

Product images:

