

Product datasheet for TP720561XL

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GADD45A (NM_001924) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human growth arrest and DNA-damage-inducible, alpha (GADD45A)

Species: Human
Expression Host: E. coli

Expression cDNA Clone

Met1-Arg165

or AA Sequence:

Tag: N-His

Predicted MW: 20.5 kDa

Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: < 0.1 EU per μg protein as determined by LAL test

Reconstitution Method: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

Stability: Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 001915

 Locus ID:
 1647

 UniProt ID:
 P24522

 Cytogenetics:
 1p31.3

Synonyms: DDIT1; GADD45





Summary: This gene is a member of a group of genes whose transcript levels are increased following

stressful growth arrest conditions and treatment with DNA-damaging agents. The protein encoded by this gene responds to environmental stresses by mediating activation of the p38/JNK pathway via MTK1/MEKK4 kinase. The DNA damage-induced transcription of this gene is mediated by both p53-dependent and -independent mechanisms. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.[provided

by RefSeq, Dec 2010]

Protein Families: Druggable Genome, Stem cell - Pluripotency

Protein Pathways: Cell cycle, MAPK signaling pathway, p53 signaling pathway

Product images:

