

## Product datasheet for **TP720618**

### **PMP2 (NM\_002677) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Human peripheral myelin protein 2 (PMP2)
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	Met1-Val132
<b>Tag:</b>	N-His
<b>Predicted MW:</b>	17.07 kDa
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>95% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
<b>Endotoxin:</b>	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
<b>Reconstitution Method:</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 μg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
<b>Storage:</b>	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
<b>Stability:</b>	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
<b>RefSeq:</b>	<u><a href="#">NP_002668</a></u>
<b>Locus ID:</b>	5375
<b>UniProt ID:</b>	<u><a href="#">P02689</a></u>
<b>RefSeq Size:</b>	3579
<b>Cytogenetics:</b>	8q21.13
<b>RefSeq ORF:</b>	396
<b>Synonyms:</b>	CMT1G; FABP8; M-FABP; MP2; P2



[View online »](#)

**Summary:**

The protein encoded by this gene localizes to myelin sheaths of the peripheral nervous system. The encoded protein can bind both the membrane layers of the sheaths and monomeric lipids, and is thought to provide stability to the sheath. A defect in this gene was shown to be a cause of dominant demyelinating CMT neuropathy. [provided by RefSeq, Jan 2017]