

## Product datasheet for **TP720634XL**

### CD46 (NM\_153826) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human CD46 molecule, complement regulatory protein (CD46), transcript variant d
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Cys35-Asp328
Tag:	C-His
Predicted MW:	33.83 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_722548</a>
Locus ID:	4179
UniProt ID:	<a href="#">P15529</a>
RefSeq Size:	3264
Cytogenetics:	1q32.2
RefSeq ORF:	1152
Synonyms:	AHUS2; MCP; MIC10; TLX; TRA2.10



[View online »](#)

**Summary:** The protein encoded by this gene is a type I membrane protein and is a regulatory part of the complement system. The encoded protein has cofactor activity for inactivation of complement components C3b and C4b by serum factor I, which protects the host cell from damage by complement. In addition, the encoded protein can act as a receptor for the Edmonston strain of measles virus, human herpesvirus-6, and type IV pili of pathogenic Neisseria. Finally, the protein encoded by this gene may be involved in the fusion of the spermatozoa with the oocyte during fertilization. Mutations at this locus have been associated with susceptibility to hemolytic uremic syndrome. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jun 2010]

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Complement and coagulation cascades