

Product datasheet for **TP720642XL**

FAM3D (NM_138805) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human family with sequence similarity 3, member D (FAM3D)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Tyr26-Phe224
Tag:	C-His
Predicted MW:	23.12 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_620160
Locus ID:	131177
UniProt ID:	Q96BQ1 , A0A0A8K9B4
RefSeq Size:	1322
Cytogenetics:	3p14.2
RefSeq ORF:	672
Synonyms:	EF7; OIT1
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)