

Product datasheet for **TP720657M**

Cathepsin L (CTSL) (NM_001912) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human cathepsin L1 (CTSL1), transcript variant 1
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Thr18-Val333
Tag:	C-His
Predicted MW:	36.88 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of 20mM HAc-NaAc, 150mM NaCl, pH 4.5.
Endotoxin:	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_001903
Locus ID:	1514
UniProt ID:	P07711 , A0A024R276
RefSeq Size:	1730
Cytogenetics:	9q21.33
RefSeq ORF:	999
Synonyms:	CATL; CTSL1; MEP



[View online »](#)

Summary:

The protein encoded by this gene is a lysosomal cysteine proteinase that plays a major role in intracellular protein catabolism. Its substrates include collagen and elastin, as well as alpha-1 protease inhibitor, a major controlling element of neutrophil elastase activity. The encoded protein has been implicated in several pathologic processes, including myofibril necrosis in myopathies and in myocardial ischemia, and in the renal tubular response to proteinuria. This protein, which is a member of the peptidase C1 family, is a dimer composed of disulfide-linked heavy and light chains, both produced from a single protein precursor. Additionally, this protein cleaves the S1 subunit of the SARS-CoV-2 spike protein, which is necessary for entry of the virus into the cell. [provided by RefSeq, Aug 2020]

Protein Families:

Druggable Genome, Protease

Protein Pathways:

Antigen processing and presentation, Lysosome