

## Product datasheet for **TP720704XL**

### MASP1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human mannan-binding lectin serine peptidase 1 (MASP1)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	His20-Arg728
Tag:	C-6His
Predicted MW:	80.68 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 $\mu$ m filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/ $\mu$ g of protein (< 1 EU/ $\mu$ g)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
Locus ID:	5648
UniProt ID:	<a href="#">P48740</a>



[View online »](#)