

Product datasheet for **TP720707XL**

POMGNT1 (NM_017739) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human protein O-linked mannanose beta1,2-N-acetylglucosaminyltransferase (POMGNT1), transcript variant 1
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Leu59-Thr660
Tag:	C-His
Predicted MW:	69.3 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_060209
Locus ID:	55624
UniProt ID:	Q8WZA1
RefSeq Size:	2753
Cytogenetics:	1p34.1
RefSeq ORF:	1980
Synonyms:	gnT-I.2; GNTI.2; GnT I.2; LGMD2O; LGMDR15; MEB; MGAT1.2; RP76
Summary:	This gene encodes a type II transmembrane protein that resides in the Golgi apparatus. It participates in O-mannosyl glycosylation and is specific for alpha linked terminal mannanose. Mutations in this gene may be associated with muscle-eye-brain disease and several congenital muscular dystrophies. Alternatively spliced transcript variants that encode different protein isoforms have been described. [provided by RefSeq, Feb 2014]



[View online »](#)

Protein Pathways: O-Mannosyl glycan biosynthesis