

Product datasheet for **TP720719XL**

UBE2B (NM_003337) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human ubiquitin-conjugating enzyme E2B (UBE2B)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Met1-Ser152
Tag:	C-His
Predicted MW:	18.3 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 μ m filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/ μ g of protein (< 1 EU/ μ g)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_003328
Locus ID:	7320
UniProt ID:	P63146
RefSeq Size:	2631
Cytogenetics:	5q31.1
RefSeq ORF:	456
Synonyms:	E2-17kDa; HHR6B; HR6B; RAD6B; UBC2



[View online »](#)

Summary: The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is required for post-replicative DNA damage repair. Its protein sequence is 100% identical to the mouse, rat, and rabbit homologs, which indicates that this enzyme is highly conserved in eukaryotic evolution. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Ubiquitin mediated proteolysis