

## Product datasheet for **TP720723**

### **PRG3 (NM\_006093) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Human proteoglycan 3 (PRG3)
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293
<b>Expression cDNA Clone or AA Sequence:</b>	Leu18-Phe225
<b>Tag:</b>	C-His
<b>Predicted MW:</b>	24.6 kDa
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>95% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
<b>Endotoxin:</b>	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
<b>Reconstitution Method:</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 μg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
<b>Storage:</b>	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
<b>Stability:</b>	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
<b>RefSeq:</b>	<u><a href="#">NP_006084</a></u>
<b>Locus ID:</b>	10394
<b>UniProt ID:</b>	<u><a href="#">Q9Y2Y8</a></u>
<b>RefSeq Size:</b>	868
<b>Cytogenetics:</b>	11q12.1
<b>RefSeq ORF:</b>	675
<b>Synonyms:</b>	MBP2; MBPH


[View online »](#)

**Summary:**

Possesses similar cytotoxic and cytostimulatory activities to PRG2/MBP. In vitro, stimulates neutrophil superoxide production and IL8 release, and histamine and leukotriene C4 release from basophils.[UniProtKB/Swiss-Prot Function]