

Product datasheet for **TP720727XL**

Kallistatin (SERPINA4) (NM_006215) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 4 (SERPINA4)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Gln21-Pro427
Tag:	C-His
Predicted MW:	47.1 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_006206
Locus ID:	5267
UniProt ID:	P29622 , A0A024R6I9
RefSeq Size:	1700
Cytogenetics:	14q32.13
RefSeq ORF:	1281
Synonyms:	KAL; kallistatin; KLST; KST; PI-4; PI4



[View online »](#)

Summary: Inhibits human amidolytic and kininogenase activities of tissue kallikrein. Inhibition is achieved by formation of an equimolar, heat- and SDS-stable complex between the inhibitor and the enzyme, and generation of a small C-terminal fragment of the inhibitor due to cleavage at the reactive site by tissue kallikrein.[UniProtKB/Swiss-Prot Function]

Protein Families: Druggable Genome, Secreted Protein