

Product datasheet for TP720737XL

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PD L2 (PDCD1LG2) (NM_025239) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human programmed cell death 1 ligand 2 (PDCD1LG2)

Species: Human
Expression Host: HEK293

Expression cDNA Clone

Leu20-Pro219

or AA Sequence:

Tag: C-His

Predicted MW: 23.59 kDa

Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)

Reconstitution Method: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

Stability: Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

819

RefSeq: NP 079515

 Locus ID:
 80380

 UniProt ID:
 Q9BQ51

 RefSeq Size:
 1518

 Cytogenetics:
 9p24.1

RefSeq ORF:

Synonyms: B7DC; bA574F11.2; Btdc; CD273; PD-L2; PDCD1L2; PDL2





PD L2 (PDCD1LG2) (NM_025239) Human Recombinant Protein - TP720737XL

Summary: Involved in the costimulatory signal, essential for T-cell proliferation and IFNG production in a

PDCD1-independent manner. Interaction with PDCD1 inhibits T-cell proliferation by blocking cell cycle progression and cytokine production (By similarity).[UniProtKB/Swiss-Prot Function]

Protein Families: Transmembrane

Protein Pathways: Cell adhesion molecules (CAMs)