

Product datasheet for **TP720741L**

SGSH (NM_000199) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human N-sulfoglucosamine sulfohydrolase (SGSH)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Arg21-Leu502
Tag:	C-His
Predicted MW:	55.72 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 μ m filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/ μ g of protein (< 1 EU/ μ g)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_000190
Locus ID:	6448
UniProt ID:	P51688
RefSeq Size:	2770
Cytogenetics:	17q25.3
RefSeq ORF:	1506
Synonyms:	HSS; MPS3A; SFMD
Summary:	This gene encodes the enzyme sulfamidase; one of several enzymes involved in the lysosomal degradation of heparan sulfate. Mutations in this gene are associated with the lysosomal storage disease mucopolysaccharidosis IIIA, also known as Sanfilippo syndrome A, which results from impaired degradation of heparan sulfate. Transcripts of varying sizes have been reported but their biological validity has not been determined. [provided by RefSeq, Jun 2017]



[View online »](#)

Protein Families: Druggable Genome

Protein Pathways: Glycosaminoglycan degradation, Lysosome, Metabolic pathways