

## Product datasheet for **TP720854L**

### MAP1LC3A (NM\_032514) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human microtubule-associated protein 1 light chain 3 alpha (MAP1LC3A), transcript variant 1
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Phe121
Tag:	C-His
Predicted MW:	15.3 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_115903</a>
Locus ID:	84557
UniProt ID:	<a href="#">Q9H492</a>
RefSeq Size:	1048
Cytogenetics:	20q11.22
RefSeq ORF:	363
Synonyms:	ATG8E; LC3; LC3A; MAP1ALC3; MAP1BLC3



[View online »](#)

**Summary:**

MAP1A and MAP1B are microtubule-associated proteins which mediate the physical interactions between microtubules and components of the cytoskeleton. MAP1A and MAP1B each consist of a heavy chain subunit and multiple light chain subunits. The protein encoded by this gene is one of the light chain subunits and can associate with either MAP1A or MAP1B. Two transcript variants encoding different isoforms have been found for this gene. The expression of variant 1 is suppressed in many tumor cell lines, suggesting that may be involved in carcinogenesis. [provided by RefSeq, Feb 2012]