

## Product datasheet for **TP720855XL**

### GRP (LGALSL) (NM\_014181) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human galectin-related protein (HSPC159)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ala2-Ser172
Tag:	C-His
Predicted MW:	20 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, pH 8.0.
Endotoxin:	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_054900</a>
Locus ID:	29094
UniProt ID:	<a href="#">Q3ZCW2</a>
RefSeq Size:	3842
Cytogenetics:	2p14
RefSeq ORF:	516
Synonyms:	GRP; HSPC159
Summary:	Does not bind lactose, and may not bind carbohydrates.[UniProtKB/Swiss-Prot Function]



[View online »](#)