

## Product datasheet for **TP720856XL**

### Parvin alpha (PARVA) (NM\_018222) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human parvin, alpha (PARVA)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Glu372
Tag:	C-His
Predicted MW:	43.2 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of 50mM Tris-HCl, 150mM NaCl, 40% Glycerol, pH 7.4.
Endotoxin:	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_060692</a>
Locus ID:	55742
UniProt ID:	<a href="#">Q9NVD7</a>
RefSeq Size:	8719
Cytogenetics:	11p15.3
RefSeq ORF:	1116
Synonyms:	CH-ILKBP; MXRA2
Summary:	This gene encodes a member of the parvin family of actin-binding proteins. Parvins are associated with focal contacts and contain calponin homology domains that bind to actin filaments. The encoded protein is part of the integrin-linked kinase signaling complex and plays a role in cell adhesion, motility and survival. [provided by RefSeq, Dec 2010]
Protein Pathways:	Focal adhesion



[View online »](#)