

Product datasheet for TP720859M

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SULT2B1 (NM 177973) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human sulfotransferase family, cytosolic, 2B, member 1

(SULT2B1), transcript variant 2

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

Met1-Glu311

Tag: C-His

Predicted MW: 36.6 kDa

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 814444

 Locus ID:
 6820

 UniProt ID:
 000204

 RefSeq Size:
 1228

Cytogenetics: 19q13.33

RefSeq ORF: 1095

Synonyms: ARCI14; HSST2



SULT2B1 (NM_177973) Human Recombinant Protein - TP720859M

Summary: Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones,

 $neurotransmitters,\,drugs,\,and\,xenobiotic\,compounds.\,These\,cytosolic\,enzymes\,are\,different$

in their tissue distributions and substrate specificities. The gene structure (number and

length of exons) is similar among family members. This gene sulfates

dehydroepiandrosterone but not 4-nitrophenol, a typical substrate for the phenol and estrogen sulfotransferase subfamilies. Two alternatively spliced variants that encode

different isoforms have been described. [provided by RefSeq, Jul 2008]

Protein Pathways: Androgen and estrogen metabolism, Sulfur metabolism