

## Product datasheet for TP720870L

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## PADI4 (NM 012387) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Human peptidyl arginine deiminase, type IV (PADI4)

Species: Human
Expression Host: E. coli

**Expression cDNA Clone** 

Met1-Pro663

or AA Sequence:

Tag: N-His

**Predicted MW:** 76.2 kDa

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

**Endotoxin:** Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

**RefSeq:** NP 036519

 Locus ID:
 23569

 UniProt ID:
 Q9UM07

 RefSeq Size:
 2265

 Cytogenetics:
 1p36.13

RefSeq ORF: 1989

Synonyms: PAD; PAD4; PAD15; PD14; PD15

Summary: This gene is a member of a gene family which encodes enzymes responsible for the

conversion of arginine residues to citrulline residues. This gene may play a role in granulocyte and macrophage development leading to inflammation and immune response. [provided by

RefSeq, Jul 2008]

