

## Product datasheet for **TP720877L**

### SAA4 (NM\_006512) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human serum amyloid A4, constitutive (SAA4)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Glu19-Thr130
Tag:	C-His
Predicted MW:	14 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 $\mu$ m filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/ $\mu$ g of protein (< 1 EU/ $\mu$ g)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_006503</a>
Locus ID:	6291
UniProt ID:	<a href="#">P35542</a>
RefSeq Size:	675
Cytogenetics:	11p15.1
RefSeq ORF:	390
Synonyms:	C-SAA; CSAA
Summary:	Major acute phase reactant. Apolipoprotein of the HDL complex.[UniProtKB/Swiss-Prot Function]



[View online »](#)