

Product datasheet for **TP720879XL**

TRIM (TRAT1) (NM_016388) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human T cell receptor associated transmembrane adaptor 1 (TRAT1)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Asn29-Asn186
Tag:	C-His
Predicted MW:	19.4 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_057472
Locus ID:	50852
UniProt ID:	Q6PIZ9
RefSeq Size:	1695
Cytogenetics:	3q13.13
RefSeq ORF:	558
Synonyms:	HSPC062; pp29/30; TCRIM; TRIM
Summary:	Stabilizes the TCR (T-cell antigen receptor)/CD3 complex at the surface of T-cells. [UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)