

Product datasheet for TP720890M

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

UBE2I (NM 194261) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human ubiquitin-conjugating enzyme E2I (UBE2I), transcript

variant 4

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

Met1-Ser158

Tag: N-GST

Predicted MW: 44.4 kDa

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 919237

 Locus ID:
 7329

 UniProt ID:
 P63279

 RefSeq Size:
 2843

Cytogenetics: 16p13.3

RefSeq ORF: 474

Synonyms: C358B7.1; P18; UBC9



UBE2I (NM_194261) Human Recombinant Protein - TP720890M

Summary: The modification of proteins with ubiquitin is an important cellular mechanism for targeting

abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. Four alternatively spliced transcript variants encoding the same

protein have been found for this gene. [provided by RefSeq, Jul 2008]

Protein Pathways: Ubiquitin mediated proteolysis