

Product datasheet for TP720900XL

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

LYPLA2 (NM_007260) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human lysophospholipase II (LYPLA2)

Species: Human
Expression Host: E. coli

Expression cDNA Clone

Met1-Val231

or AA Sequence:

Tag: C-His

Predicted MW: 25.8 kDa

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 009191

Locus ID: 11313

UniProt ID: 095372, A0A140V|C9

RefSeq Size: 1648
Cytogenetics: 1p36.11
RefSeq ORF: 693

Synonyms: APT-2; APT2; DJ886K2.4

Summary: Lysophospholipases are enzymes that act on biological membranes to regulate the

multifunctional lysophospholipids. There are alternatively spliced transcript variants

described for this gene but the full length nature is not known yet. [provided by RefSeq, Jul

2008]

Protein Pathways: Glycerophospholipid metabolism

