

Product datasheet for TP720913

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

P4HB (NM_000918) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human prolyl 4-hydroxylase, beta polypeptide (P4HB)

Species: Human
Expression Host: E. coli

Expression cDNA Clone

Asp18-Lys505

or AA Sequence:

Tag: C-His

Predicted MW: 55.9 kDa

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 000909

Locus ID: 5034

UniProt ID: <u>P07237</u>, <u>A0A024R8S5</u>

RefSeq Size: 2596 Cytogenetics: 17q25.3 RefSeq ORF: 1524

Synonyms: CLCRP1; DSI; ERBA2L; GIT; P4Hbeta; PDI; PDIA1; PHDB; PO4DB; PO4HB; PROHB



ORIGENE

Summary: This gene encodes the beta subunit of prolyl 4-hydroxylase, a highly abundant

multifunctional enzyme that belongs to the protein disulfide isomerase family. When present as a tetramer consisting of two alpha and two beta subunits, this enzyme is involved in hydroxylation of prolyl residues in preprocollagen. This enzyme is also a disulfide isomerase

containing two thioredoxin domains that catalyze the formation, breakage and rearrangement of disulfide bonds. Other known functions include its ability to act as a

chaperone that inhibits aggregation of misfolded proteins in a concentration-dependent manner, its ability to bind thyroid hormone, its role in both the influx and efflux of S-nitrosothiol-bound nitric oxide, and its function as a subunit of the microsomal triglyceride

transfer protein complex. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome