

Product datasheet for **TP720928XL**

SULT4A1 (NM_014351) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human sulfotransferase family 4A, member 1 (SULT4A1)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Leu284
Tag:	Tag Free
Predicted MW:	33 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Endotoxin:	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_055166
Locus ID:	25830
UniProt ID:	Q9BR01
RefSeq Size:	2465
Cytogenetics:	22q13.31
RefSeq ORF:	849
Synonyms:	BR-STL-1; BRSTL1; DJ388M5.3; hBR-STL-1; NST; SULTX3
Summary:	This gene encodes a member of the sulfotransferase family. The encoded protein is a brain-specific sulfotransferase believed to be involved in the metabolism of neurotransmitters. Polymorphisms in this gene may be associated with susceptibility to schizophrenia. [provided by RefSeq, Jul 2008]



[View online »](#)