

Product datasheet for **TP720929XL**

Hsp40 (DNAJB1) (NM_006145) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human DnaJ (Hsp40) homolog, subfamily B, member 1 (DNAJB1)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Gly2-Ile340
Tag:	C-His
Predicted MW:	39.1 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_006136
Locus ID:	3337
UniProt ID:	P25685 , Q6FHS4
RefSeq Size:	2419
Cytogenetics:	19p13.12
RefSeq ORF:	1020
Synonyms:	Hdj1; Hsp40; HSPF1; RSPH16B; Sis1



[View online »](#)

Summary:

This gene encodes a member of the DnaJ or Hsp40 (heat shock protein 40 kD) family of proteins. DNAJ family members are characterized by a highly conserved amino acid stretch called the 'J-domain' and function as one of the two major classes of molecular chaperones involved in a wide range of cellular events, such as protein folding and oligomeric protein complex assembly. The encoded protein is a molecular chaperone that stimulates the ATPase activity of Hsp70 heat-shock proteins in order to promote protein folding and prevent misfolded protein aggregation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015]