

## Product datasheet for TP720932XL

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## YKT6 (NM 006555) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Human YKT6 v-SNARE homolog (S. cerevisiae) (YKT6)

Species: Human
Expression Host: E. coli

**Expression cDNA Clone** 

Met1-Met198

or AA Sequence:

Tag: N-His

**Predicted MW:** 24.5 kDa

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

**RefSeq:** NP 006546

**Locus ID:** 10652

UniProt ID: <u>015498</u>, <u>A4D2|0</u>

RefSeq Size: 2783
Cytogenetics: 7p13
RefSeq ORF: 594

Summary: This gene product is one of the SNARE recognition molecules implicated in vesicular transport

between secretory compartments. It is a membrane associated, isoprenylated protein that functions at the endoplasmic reticulum-Golgi transport step. This protein is highly conserved from yeast to human and can functionally complement the loss of the yeast homolog in the

yeast secretory pathway. [provided by RefSeq, Jul 2008]

**Protein Pathways:** SNARE interactions in vesicular transport

