

Product datasheet for **TP720961XL**

Cxcl16 (NM_023158) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse chemokine (C-X-C motif) ligand 16 (Cxcl16)
Species:	Mouse
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Met1-Lys218
Tag:	tag free
Predicted MW:	20.1 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_075647
Locus ID:	66102
UniProt ID:	Q8BSU2 , A2CFE9
RefSeq Size:	2448
Cytogenetics:	11 B3
RefSeq ORF:	741
Synonyms:	0910001K24Rik; AV290116; b2b498Clo; BB024863; CXCL16v1; CXCL16v2; SR-PSOX; Zmynd15



[View online »](#)

Summary:

Induces a strong chemotactic response. Induces calcium mobilization. Binds to CXCR6/Bonzo. Also acts as a scavenger receptor on macrophages, which specifically binds to OxLDL (oxidized low density lipoprotein), suggesting that it may be involved in pathophysiology such as atherogenesis.[UniProtKB/Swiss-Prot Function]