

Product datasheet for TP720972XL

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

UBE2Z (NM 023079) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Human ubiquitin-conjugating enzyme E2Z (UBE2Z)

Species: Human
Expression Host: E. coli

Expression cDNA Clone

Met1-Val246

or AA Sequence:

Tag: N-His

Predicted MW: 30.2 kDa

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 075567

 Locus ID:
 65264

 UniProt ID:
 Q9H832

 RefSeq Size:
 3101

Cytogenetics: 17q21.32

RefSeq ORF: 738

Synonyms: HOYS7; USE1

Summary: This gene encodes an enzyme which ubiquitinates proteins which participate in signaling

pathways and apoptosis. [provided by RefSeq, Feb 2012]

Protein Pathways: Acute myeloid leukemia, Pathways in cancer, Ubiquitin mediated proteolysis

