

Product datasheet for **TP720981**

PIN4 (NM_006223) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human protein (peptidylprolyl cis/trans isomerase) NIMA-interacting, 4 (parvulin) (PIN4), nuclear gene encoding mitochondrial protein, transcript variant 1
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Lys156
Tag:	N-His
Predicted MW:	18.8 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	Endotoxin level is < 0.1 ng/µg of protein (< 1 EU/µg)
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_006214
Locus ID:	5303
UniProt ID:	Q9Y237
RefSeq Size:	1331
Cytogenetics:	Xq13.1
RefSeq ORF:	468
Synonyms:	EPVH; PAR14; PAR17



[View online »](#)

Summary:

This gene encodes a member of the parvulin subfamily of the peptidyl-prolyl cis/trans isomerase protein family. The encoded protein catalyzes the isomerization of peptidylprolyl bonds, and may play a role in the cell cycle, chromatin remodeling, and/or ribosome biogenesis. The encoded protein may play an additional role in the mitochondria. [provided by RefSeq, Dec 2009]