

## Product datasheet for **TP720983L**

### Syntenin (SDCBP) (NM\_001007070) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human syndecan binding protein (syntenin) (SDCBP), transcript variant 5
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ser2-Val298
Tag:	C-His
Predicted MW:	33 kDa
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 40% glycerol, PH8.0.
Endotoxin:	Endotoxin level is < 0.1 ng/μg of protein (< 1 EU/μg)
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001007071</a>
Locus ID:	6386
UniProt ID:	<a href="#">O00560</a>
RefSeq Size:	2167
Cytogenetics:	8q12.1
RefSeq ORF:	891
Synonyms:	MDA-9; MDA9; ST1; SYCL; TACIP18



[View online »](#)

**Summary:**

The protein encoded by this gene was initially identified as a molecule linking syndecan-mediated signaling to the cytoskeleton. The syntenin protein contains tandemly repeated PDZ domains that bind the cytoplasmic, C-terminal domains of a variety of transmembrane proteins. This protein may also affect cytoskeletal-membrane organization, cell adhesion, protein trafficking, and the activation of transcription factors. The protein is primarily localized to membrane-associated adherens junctions and focal adhesions but is also found at the endoplasmic reticulum and nucleus. Alternative splicing results in multiple transcript variants encoding different isoforms. Related pseudogenes have been identified on multiple chromosomes. [provided by RefSeq, Jan 2017]

**Protein Families:**

Druggable Genome, Transmembrane